Are you prepared for safe emergency evacuation?

If you are not sure, then we can help.



THE WORLD'S NO.1 STAIRWAY EVACUATION CHAIR SINCE 1982

What makes us better than the rest?



- World's #1 EVACuation chair
- Single person operation*
- Lightweight, quick and easy to use
- 500lb weight capacity**
- Versatile Can also be used along corridors, hallways and flat surfaces
- Folds away Ready to use in seconds
- Wall mounted with dust cover
- 🕂 Longest warranty Lifetime
- ✤ Over 35 years experience
- ✤ Never a product failure or recall
- Successful evacuations during 9/11
- FDA Registered
- Exceeds fire safety requirements and emergency response guidelines established to support the Americans with Disabilities Act

*Additional assistance and practice is recommended for heavy passengers **500lb weight capacity applies to Evac+Chair 500H

Our pledge

OUR PRICE GUARANTEE

Our price promise guarantees if you find a similar product for a better price or the same product cheaper elsewhere, we will meet or beat that price.

LIFE TIME SUPPORT

We are so confident our products are the best on the market that we offer a lifetime manufacturer's warranty with every chair sold*.

CUPPORA * LIFE TIME * SUPPOR

NEVER COMPROMISE ON QUALITY

You have complete peace of mind knowing your chair is guaranteed to be authentic and the best possible quality. Our certified products are rigorously tested to exceed all fire, safety and emergency response guidelines.

100% CUSTOMER SATISFACTION

We aim to deliver superior service, every time, with 100% satisfaction guaranteed. As the world's number 1 and the original Evacuation Chair company, we have more satisfied customers than anyone else.





*Lifetime warranty excludes Evac+Chair 700H and Evac+Power 800.

What are your current emergency evacuation plans?

How will everyone escape?

- Most people at some time during their lives will have a disability, either temporary or permanent, that will limit their ability to evacuate on their own.
- FACT: More than one in five Americans have some type of disability.

Waiting for first responders

- Waiting with the person with the impairment for first responders should be the last choice in an emergency.
- SECONDS COUNT: There is no way of predicting how long an area of safe refuge will remain safe.

Can the elevators be used?

- Elevators stop working when an alarm is activated.
- Restrictions are imposed on the use of elevators during some types of building emergencies.

Mobility impairments

 Anyone with a mobility issue will be unable to evacuate safely unless proper evacuation equipment is provided.

Other impairments and multiple impairments

- Anyone at anytime might need the assistance of evacuation products in an emergency.
- Broken bones, illness, trauma, pregnancy, heart or lung condition, panic attack, arthritis, neurological diseases and rheumatism can all be an impairment to evacuation.
- Reduction in overall ability is also experienced by people as they age. People of extreme size or weight (Bariatric) often need assistance as well.



Evac+Chair® 300H

Evac+Chair[®] 500H

Evac+Chair[®] 600H

PRODUCT OVERVIEW

The EVAC+CHAIR[®], model 300H (MK4) is a device which permits the safe egress of mobility impaired persons from a building in the event of an emergency or evacuation. The standard model is our lightest weighing stair chair at only 19lbs. The 300H features a 400lb. weight capacity and a blue power coated textured finish with a contrasting fire retardant vellow hammock.

DUAL

POSITION

SEAT

KEY FEATURES

- Adjustable patient head restraint
 Friction tracks
- Ergonomic adjustable handle
- Head, chest and leg restraints
- Lockable rear wheels
- ACCESSORIES INCLUDED
- LIFETIME WARRANTY
- Photoluminescent signage
- Wall mounting brackets
- Instructional DVD
- Vinvl dust cover
- User guide

SPECIFICATION

41" (folded) Height

- Width 20" (folded)
- 8" (folded) Depth
- 19lbs Weight
- Weight Capacity 400lbs

PRODUCT OVERVIEW

The EVAC+CHAIR[®], 500H is a device which permits the safe egress of mobility impaired persons from a building in the event of an emergency or evacuation. A one person operation enhanced EVAC+CHAIR[®] with an increased passenger payload capacity of 500lbs. and for use on steep stairs (45°).

KEY FEATURES

Adjustable patient head restraint
 Friction tracks

500lb

WEIGHT CAPACITY

- Ergonomic adjustable handle
- Fold away grip handle

ACCESSORIES INCLUDED

LIFETIME WARRANTY

- Photoluminescent signage
- Wall mounting brackets
- User guide

SPECIFICATION Height 45" (folded)

Width 20" (folded) Depth 9" (folded) Weight 29lbs Weight Capacity 500lbs

PRODUCT OVERVIEW

The EVAC+CHAIR[®]. 600H is a device which permits the safe egress of mobility impaired persons from a building in the event of an emergency or evacuation. The EVAC+CHAIR® 600H incorporates 2 carry-handles

that allows for two person operation for locations with difficult access, including upwards evacuation. Can also be used in single person operation.

KEY FEATURES

- Adjustable patient head restraint
 Friction tracks
- Ergonomic adjustable handle
- Head, chest and leg restraints
- Lockable rear wheels
- Fold away front and rear grab handles

ACCESSORIES INCLUDED

- LIFETIME WARRANTY
- Photoluminescent signage User guide
- Wall mounting brackets
- Instructional DVD



Vinyl dust cover

Reflective decals

Heavy duty wheels

Dual position seat

EVACOCHAIR

SPECIFICATION

Height	47" (folded)
Width	23" (folded)
Depth	12" (folded)
Weight	24lbs
Weight Capacity	400lbs

- Vinyl dust cover

- Reflective decals
- Heavy duty wheels

EVACOCHAID

- Dual position seat
- 4 patient safety straps

I ockable rear wheels

- Heavy duty wheels Dual position seat



EVACOCHAI

Reflective decals

Evac+Chair® 700H

Evac+Power 800

Carrylite Transit Chair

PRODUCT OVERVIEW

The EVAC+CHAIR[®], IBEX Tran-Seat safely ascends and descends stairs without carrying or heavy lifting. The IBEX 700H requires 2 person operation and reduces the physical demands on the operators and reduces the risk of muscle strain, hernia and back injury.

Adjustable brake

KEY FEATURES

- Telescopic extendable handle
- Stainless steel frame construction Footrest
- Detachable seat and backrest
- Adjustable patient belts
- Front guide handles

ACCESSORIES INCLUDED

- LIFETIME WARRANTY
- Patient Restraints
- Integrated handles with no pinch points

SCENDS

- Wall mounting brackets
- Instructional DVD
- Vinyl dust cover

PRODUCT OVERVIEW

This battery-powered, electric evacuation chair can be used in emergency evacuations, power outages, and daily patient transport to safely and easily move a passenger up or down multiple flights of stairs without the need for carrying or lifting.

SCENDS

AND

DESCENDS

KEY FEATURES

- Quick charge lithium ion battery pack
- High strength aluminum construction
- Adjustable stair tracks
- LED safety light
- Quick release four point harness

ACCESSORIES INCLUDED

MANUFACTURERS WARRANTY

45" (folded)

19.5" (folded)

10.6" (folded)

66lbs

350lbs

- Battery Charger (US/EU/UK)
- Rechargeable Battery
- Photoluminecent Wall Sign
- User Guide

SPECIFICATION

Weight Capacity

Height

Width

Depth

Weight

PRODUCT OVERVIEW

Designed for easy storage the EVAC+CHAIR[®], CarryLite Transit Chair folds neatly in half and can be stored in the smallest of areas. The chair unfolds in one swift movement for immediate use.

KEY FEATURES

- Lightweight aluminum construction
- Integrated handles (no pinch points)
- Unique locking system
- Adjustable patient straps
- Large wheels
- Footrest

ACCESSORIES INCLUDED

- LIFETIME WARRANTY
- Passenger safety belt



SPECIFICATIONHeight43" (folded)Width18" (folded)Depth14" (folded)Weight32lbsWeight Capacity350lbs

Evacuation Sleds ENABLING THE MOBILITY IMPAIRED TO 'SLIDE TO SAFETY'

ResQme[®]



Product Overview

Quick and easy to prepare for evacuation, allows ONE person to safely transport another in emergencies. Easily moves over ANY surface, across floors and down stairs. It is the ONLY Mat that easily adjusts to the passenger's size (one size fits ALL). Soft-absorbing head pillow reduces head movement going down stairs. Long towing straps at the head and feet for easy maneuvering.

Key Features

- Pulling Straps at both ends
- Head pillow for comfort and protection
- Will slide on carpet, tiles vinyl and concrete
- Can be pulled by one person
- Strong 4" wide velcro straps to secure the person being rescued
- Wrist/Arm strap to stop hands and arms from flaving out of control
- Can be used for children and adults up to 6'6" tall
- Base board to protect the person from hitting door frames etc., when going around corners

ResQmat[®]



Product Overview

The ResQmat[®] is the perfect solution for the swift emergency evacuation of persons who are mobility impaired. UNROLL STRAP and GO. The ResQmat is used ALL over the world in Hospitals, Nursing Homes, Retirement Homes, Assisted Living, and Ambulance Companies....ANYWHERE!

Key Features

- Quick and Easy to prepare for evacuation
- Allows ONE person to safely transport another in emergencies
- Easily moves over ANY surface across floors and down stairs
- ONLY Mat that easily adjusts to the passenger's size. (one size fits ALL)
- Soft-absorbing head pillow reduces head movement going down stairs
- Long towing straps at the head and feet for easy maneuvering

ResQpod[®]



The ResQpod is an evacuation board designed to help evacuate non-ambulant people or people who may not be able to walk down stairs in an emergency. This high-performance evacuation board can quickly and safely get children and adults to safety. The ResQpod can be used in office buildings, theatres, libraries, hospitals, schools, and many other facilities. Additionally, this evacuation board comes with a lifetime warranty!

Key Features

- Long pulling straps located at both head and feet to allow for multiple users
- Secure body protection which completely encloses the upper body
- Heavy duty hook and loop for increased passenger safety
- Impact absorbing pillow for added head protection when transporting down stairs or across rough terrain
- Glideable plastic board corrugated plastic board allows for excellent heat insulation preventing injury, plus removing the necessity to unroll every few months to prevent the plastic keeping its shape
- Adjustable straps curl the ResQPod around the feet, allowing it to fit all sizes

Rescue Sheet



Product Overview

The Rescue Sheet is widely accepted as a vital aid in the rapid evacuation of nonambulatory patients during emergency situations. The Rescue Sheet is located permanently and unobtrusively under the mattress and does not interfere with bed making. The patient is rescued with their bed mattress to offer protection under the patient. especially when going down stairs.

Please Note: Also Available for Vinyl & Carpet

Key Features

- Permanently in place between the mattress and bed frame
- Fits all bed models and mattresses
- Ready for immediate use
- Long straps located at the head and feet for easy maneuvering
- Adjustable wrap around straps for safely cocooning the patient
- Heavy duty construction for maximum efficiency
- Designed to enable one or more persons to evacuate a person in their mattress during an emergency



Accessories



Comfy Seat (300H, 500H, 600H)

A cushioned seat which provides additional comfort and support by clipping to the seat frame of the EVAC+CHAIR[®] and can be purchased as a retrofit.

Anti-Theft Alarm Device

The Anti-Theft Alarm device batteryoperated siren which gives out an audible alarm when the EVAC+CHAIR[®] is detached from its housing.





iEvac®Smoke

The iEvac® protects against byproducts of combustion including carbon monoxide and hydrogen cyanide, as well as other chemical and biological agents.

Secure Steel Cabinet

A 20 gauge galvanized steel storage cabinet. Has a wide manual door, with nonlocking paddle handle, for easy access and quick removal (Optional locking handle).





Evac+Chair® Stand

A self standing support suitable for mounting the Evac+Chair avoiding the need for wall mounting brackets. Internal use only.

Patient Carrier Seat

A portable transfer device for lifting incapacitated persons from a wheelchair into or out of an EVAC+CHAIR[®], where self-transfer is difficult.



How does legislation affect your business?

- NFPA 101, 2015
- The Rehabilitation Act of 1973
- Disability Discrimination & Reasonable Accommodation
- Americans with Disabilities Act of 1990
- Part 1630 Regulations to implement the equal employment provisions of the Americans with disabilities Act
- Civil Rights Act of 1991
- United States Access Board Required Compliance with the IBC
- International Building Code (IBC)-2000
- American National Standards for Building and Facilities
- Title III
- The Disabled Persons Act
- Unruh Act (California Civil Code Section 51)
- New York City Human Rights Law
- Maryland House Bill 1061
- California Assembly Bill 3098



EVAC^C**HAIR**[®]

THE WORLD'S NO.1 STAIRWAY EVACUATION CHAIR SINCE 1982



EVAC+CHAIR® AND ITS SILHOUETTE SYMBOL ARE REGISTERED TRADEMARKS. COPYRIGHT BY EVAC+CHAIR® NORTH AMERICA LLC. ALL RIGHTS RESERVED. For further information please contact:

Evac+Chair North America

3000 Marcus Ave, Suite #3E6, Lake Success, NY, 11042-1012

> T 516-502-4240 F 516-327-8220 W www.evac-chair.com E sales@evac-chair.com